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### GB Product information for surgical instruments

#### Contents:

1. Product description
2. Indications of surgical instruments
3. Contraindications of surgical instruments
4. Notes on processing
5. Recommendations for use

#### 1. Product description

Dear customer,

The surgical instruments with the prefix RF at the beginning of the order number are made of stainless steel. The instruments with the prefix H at the beginning of the order number are made of tungsten carbide. The diamond instruments are made of stainless steel and coated galvanically with diamond particles.

Please read the instructions for use carefully before using the surgical instruments.

The surgical instruments were developed for use by dentists trained and experienced in surgery and may only be used by these. The users are obliged to check the material for its suitability for the intended use at their own responsibility. This is particularly necessary if this concerns a use not mentioned in the product information.

The manufacturer shall not be liable for any damages due to non-compliance with these provisions in the instructions for use or which are due to system-incompatible use.

#### 2. Indications of surgical instruments

The surgical instruments were developed for the following areas of use:

**2.1. TissueTrimmer:** is used in place of an electrotome (electrosurgery) or scalpel for cutting the gingiva. Turbine: >300.000 min<sup>-1</sup>

**Areas of use:** gingival contouring, widening the sulcus for impressions, exposing deep cervical cavities, papillectomy, recovering intraosseous implants as well as impacted teeth, removing any type of granulation tissue.

**2.2. Allport Bone Cutter:** the round cutters are used for apicectomy, work on the cortical bone, pre-drilling in bone, harvesting of bone as well as for the window technique.

**2.3. H151:** special cutter for drilling into the depth of the bone, for smoothing the bone surface.

**2.4. H161 Lindemann:** cutter for pre-drilling in bone and smoothing of bone.

**2.5. H162 Lindemann:** for separating and cutting in bone, separation of teeth and roots.

**2.6. H162A:** the design of the tip allows vertical drilling into the bone and the teeth. For separating bone, teeth and roots.

**2.7. H162S:** paddle-shaped tip for vertical drilling into bone and teeth. For separating bone, teeth and roots.

**2.8. H166 Lindemann:** for separating and cutting in bone, separation of teeth and roots.

**2.9. H166A:** the design of the tip allows vertical drilling into the bone and the teeth. For separating bone, teeth and roots.

**2.10. H167 Lindemann:** for separating and cutting in bone, separation of teeth and roots.

**2.11. H254:** for slitting of bone, through cortical and cancellous bone, for blade implants or for exposing mandibles for bone splitting.

**2.12. H254A:** the design of the tip allows vertical drilling into the bone and the teeth. For separating bone, teeth and roots. Also suitable for exposing mandibles for bone splitting.

**2.13. H255A:** the design of the tip allows vertical drilling into the bone and the teeth. For separating bone, teeth and roots. Also suitable for exposing mandibles for bone splitting.

**2.14. H267:** for separating and cutting in bone, separation of teeth and roots.

**2.15. H269:** for separating and cutting in bone, separation of teeth and roots.

**2.16. RF161 Lindemann:** for separating and cutting in bone.

**2.17. RF162 Lindemann:** for separating and cutting in bone, separation of teeth and roots

**2.18. RF163 Lindemann:** for separating and cutting in bone, separation of teeth and roots.

**2.19. Triangular drills:** these are used for pre-drilling before the pilot bur and to determine the implant axis. The depth markings on the triangular drill indicate the drilling depth (illustration).

**2.20. Trepan burs:** the hollow burs are used for harvesting autologous bone material and for explantation. Markings on the implant drills indicate the drilling depth (illustrations).

**2.21. RF210 pilot burs and RF210IK pilot burs internally cooled:** pre-drilling for implants depending on the size of the implant. After pre-drilling, it is essential to use the implant drill of the respective implant manufacturer. The markings on the implant drills indicate the drilling depth (illustrations).

**2.22. RF81IK and RF167IK - internally cooled instruments:** the RF81IK round burs are used for apicectomy, work on the cortical bone, pre-drilling in bone, harvesting of bone as well as for the window technique. The RF167IK is used for cutting and drilling in bone and for the smoothing of bone.

**2.23. 361 ChirDisc diamond disc:** bone separation in the mandible according to Prof. Khouri for bone distraction.

**2.24. CDS1 ChirDia:** separation of teeth, resecting of root remnants and pieces of bone.

**2.25. CDS2 ChirDia:** separation of teeth, suited in particular for removing deep root remnants and pieces of bone.

**2.26. 801 round diamond burs:** for the window technique in implantology, bone grinding, apicectomy, gingiva contouring for pontics.

#### 3. Contraindications of surgical instruments

The surgical instruments may under no circumstances be used outside the areas listed under Section 2. In particular, the specified speeds may not be exceeded. The instruments may not be used without sufficient coolants, exception TissueTrimmer.

#### 4. Notes on processing

The surgical instruments are supplied non-sterile and must be cleaned, disinfected and sterilised prior to initial use. Place the instruments into a cleaning and disinfection bath immediately after use. This is the only way prevent drying in of the tissue substances. Then clean, check and disinfect and sterilise as follows: use an approved aldehyde-free disinfectant VAH (D) or FDA (US) or with CE (EU) certificate. Use an agent with corrosion protection for the instruments made of stainless steel.

**Attention: Never use alcohol-based disinfectants!** Possibly with additional ultrasonic cleaning (only in rack to avoid damage), observe the concentrations according to the instructions of the manufacturer of the disinfection bath. All surgical instruments can also be cleaned in a thermal disinfectant (bur block 4063).

Check the instruments: damaged, bent and non-concentric instruments are to be sorted out immediately and no longer to be used. Then seal and sterilise in an autoclave (134°C /273°F; 2.1 bar; 5 min) and reuse.

#### Special attention for internally cooled instruments:

Use branded single-use cannulas Luer Lock No. 16.141 size 0.4x40mm for cleaning the water-conducting inner channels of the IK instruments. Then blow through with water and air using a pressure gun.

**Germany:** These instruments are classified under the rating “Critical B”. The guidelines of the Robert Koch Institute are to be observed: “Infection prevention in dentistry – Hygiene requirements” and “Hygiene requirements in the processing of medical devices”. Other than that, the legal and hygiene regulations in the respective countries shall apply!

Detailed processing recommendations according to DIN EN ISO 17664 ☛ www.nti.de

#### 5. Recommendations for use

When used in bone, the instruments may only be used with physiological saline solution (sodium chloride) 0.9%. If drilling is to be performed to greater depths, it may prove necessary to operate with additional cooling.

Only surgical contra-angle or handpiece attached to a surgical drive may be used when working in bone. The TissueTrimmer is an exception, see 5.1. Use at unsuitable angles, insufficient straight access and/or excessive force is to be avoided (risk of necroses, increased risk of fracture).

The working speeds are given on the packaging and in the catalogue and must be complied with.

**5.1. TissueTrimmer** (No. 2.1) is always used with the turbine at a speed of 300,000–500,000 min<sup>-1</sup> (a lower speed can lead to injuries and is therefore to be avoided). The cooling water is to be turned off, please check before use. The tissue surface must be dry, otherwise the TissueTrimmer cannot work effectively. The high speed of the turbine converts kinetic energy (motion energy) into heat on contact with the tissue. The working part is not heated in the process as the material is heat-resistant. The heat is immediately dissipated into the connective tissue which leads to staunching bleeding (haemostasis). This allows extremely caring and accurate working. When widening the sulcus for impressions, prepare 2-3 times slower around the tooth than when preparing with a diamond instrument. The instrument must always be used in one direction. Never work backwards and forwards. Such a procedure would prevent coagulation.

**5.2.** With the exception of the speeds, no special instructions are to be adhered to for the instruments with the numbers 2.2 to 2.19 made of tungsten carbide or stainless steel.

**5.3. Nr. 2.20. Trepan burs:** the trepan burs must first be used with anti-clockwise rotation until a groove of approx. 1.0 to 1.5 mm depth has been achieved in the bone. Then one can continue with clockwise rotation in the groove. It is essential to comply with the maximum drilling depth (determined via X-ray) to avoid injuries to the nerve pathways. In the ascending branch of the mandible, an angle of the 15°-20° to the sagittal plane is to be maintained for the trepan bur.

**5.4. Nr. 2.21.** With the exception of the speeds, no special instructions are to be adhered to for the pilot burs and the internally cooled pilot burs.

**5.5. Nr. 2.22.** With the exception of the speeds, no special instructions are to be adhered to for the internally cooled instruments RF81IK and RF167IK.

**5.6. No. 2.23. ChirDisc diamond disc:** compliance with the speeds prevents forming of a necrosis. The disc should only be used with a disc guard.

**5.7. No. 2.24. - 2.26.** With the exception of the speeds, no special instructions are to be adhered to for the diamond instruments.

**The general instructions for use for rotary dental and surgical instruments (B14-4) are to be observed! ☛ www.nti.de**

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© NTI-Kahla GmbH Rotary Dental Instruments