Slidematic Facebow
Instruction Manual

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Arco Facial Slidematic
Manual de Instrucciones

Table of Contents

I. Rationale for Development	2
II. Benefits and Features	2
III. Slidematic Facebow Procedure	4
IV. Mounting Procedure	6
V. Slidematic Accessories	8
VI. Reordering	8
VII. Warranty	8

Please check online at
www.whipmix.com
for the most current
instructions and parts list.

I. Rationale for Development

Based on in-depth analysis of the facebows on the market engineers in conjunction with several leading clinicians and educators prepared a list of criteria for the development of this new facebow. The Slidematic Facebow was designed to satisfy the following requirements:

- An accurate earbow that moves straight horizontally for different facial widths, unlike the caliper type that moves off the vertical axis.
- A minimal number of parts and finger screws to permit a facebow registration to be taken in less time.
- The ability to send the bitefork assembly (without the measuring bow) to the laboratory, preventing damage to the bow in transit.
- Ease of mounting the maxillary cast by not having to transfer the measuring bow to the articulator.
- The capability of taking multiple facebow records at the lowest possible cost through the use of only one measuring bow and different bitefork assembly for each case.

II. Benefits and Features


Facebow Transfer

The facebow transfer procedure establishes the relationship of the maxillary dentition to the horizontal reference plane so that the maxillary cast may be mounted on the articulator in the correct anatomical position. The Slidematic Facebow provides a fast, easy and extremely accurate means of transferring the proper anatomical relationship to the articulator.

Benefits and Features

- The precision manufacturing of the unique “speed-slide” gear mechanism makes it very easy to assemble the bow on the patient.
- The right and left arms of the measuring bow are geared to precise equidistant movement from the center of the bow.
- The scale on the measuring bow represents **the patient’s intercondylar distance** (not the interbow distance) for ease in setting articulators having this adjustment.
- The Slidematic Facebow can be used on all Denar® articulators. It also adapts for use with Whip Mix®, and Hanau™ articulators.
- The bow uses the external auditory meatus reference point for determining the arbitrary hinge axis location.
- The built-in reference pointer aligns the bow with the horizontal reference plane.
- The built-in sight can be used to view the anterior reference point.
- The measuring bow should not be mounted on the articulator during the transfer procedure. The bow, when detached from the vertical shaft, can be used again immediately with the additional transfer jigs (vertical shaft, bitefork assembly and articulator index). In this manner, the maxillary casts may be mounted at any time.

- In one step the bitefork assembly is secured to the articulator index and maxillary cast can be mounted to the articulator.
- All finger screws are easily accessible from the front.
- The procedures can be readily delegated to auxiliary personnel.

For a further breakdown, refer to Denar® Slidematic Facebow Parts List #8815-1.

Sterilize bitefork before each use.

III. Slidematic Facebow Procedure

Mark the anterior reference point on the patient's right side using the Reference Plane Locator and Marker. The point is 43 millimeters above the incisal edge of the right central or lateral incisors (see Figure 1). On an endentulous patient, measure up from the lower border of the upper lip when it is in repose.

Cover the metal with two thicknesses of base-plate wax softened in warm water (approximately 135°F or 55° C). With the bitefork arm to the patient's right, place the fork in the mouth, aligning the patient's midline with the index notch, so that it is parallel with the patient's coronal and horizontal planes (see Figure 2). Be certain to obtain a light indexing of the patient's maxillary arch and then ask the patient to hold the bitefork in place.

Attach the vertical shaft to the measuring bow with the clamp marked #2 on the patient's right and tighten the finger screw (see Figure 3). It is necessary to tighten this finger screw to secure the vertical shaft to the measuring bow and also to avoid movement. This same movement will occur after inserting the vertical shaft in the articulator index as shown in Figure 12. Be sure to tighten the finger screw.

Loosen the finger screws on the clamps marked #1 and #2 on the vertical shaft. With your thumb, loosen the center wheel on the top of the measuring bow and slide the bow open to accommodate the width of the patient's face. Assemble the facebow on the patient by sliding the bitefork arm through the hole in clamp #2 as the measuring bow's earpieces are placed in the patient's auditory meatus (see Figure 4). Tighten the center wheel on the measuring bow and loosen the finger screw on the anterior reference pointer. Raise or lower the bow so that the pointer or sight aligns precisely with the anterior reference point (see figure 5) and tighten clamp #1, then clamp #2 (see Figures 6 and 7). When tightening clamps #1 and #2, care must be taken

Figure 1

Figure 2


Figure 3

Figure 4

Figure 5

not to displace the bow to either side by having the vertical shaft rest on the fingers as shown in **Figure 6**. The patient's inter-condylar distance is the measurement indicated on the scale (see **Figure 8**). Record this measurement.

Loosen the finger screw on the measuring bow, slide the bow open, and remove the entire facebow from the patient.

Detach the measuring bow from the transfer jig by loosening the finger screw (see **Figure 9**). Having completed the procedures involving the patient, the bitefork assembly (see **Figure 10**) may be labeled with the patient's name and set aside while the measuring bow portion can be used with an additional bitefork assembly for the next patient.

Note: The metal bitefork, vertical shaft and earpieces can be sterilized in an autoclave EXCEPT for the black finger screws on the #1 and #2 clamps. Remove the finger screws (and spacers) before autoclaving or use cold sterilization.

Sterilize bitefork before each use.

Figure 6

Figure 7


Figure 8

Figure 9

Figure 10

IV. Mounting Procedure

A benefit of using the Denar® Slidematic Facebow is that multiple transfer jigs may be used with only one measuring bow. In some cases the mounting of the maxillary cast can be delegated to the laboratory, involving no loss of accuracy and no period of time without facebow transfer capability in the dental office. The laboratory can attach an articulator index to its own Denar® articulator and mount the maxillary cast using only the bitefork assembly from the dental office. Each articulator index positions the bitefork assembly on any Denar® articulator so that the relationship with the condyles recorded on the patient is accurately reproduced on the articulator.

Denar Articulator Index

Replace the incisal table on the articulator with the articulator index (see Figure 11). With the numbers on clamps #1 and #2 in the upright position, secure the vertical shaft of the bitefork assembly in the hole of the articulator index. Tighten the finger screw on the front edge of the index (see Figure 12). The incisal pin should set on top of the slide insert on the articulator index (see Figure 13).

Hanau Articulator Index

Remove the mounting plate from the lower member. Using the lower member thumb screw attach the Hanau index with the transfer assembly hole to the front of the articulator. With the numbers on clamps #1 and #2 in the upright position, secure the vertical shaft of the bitefork assembly in the hole of the Hanau index. Tighten the finger screw on the front edge of the index (see Figure 14).

Whip Mix Screw-type Index

Remove the mounting plate from the lower member. Using the lower member thumb screw attach the Whip Mix screw-type index with the transfer assembly hole to the front of the articulator. With the numbers on clamps #1 and #2 in the upright position, secure the vertical shaft of the bitefork assembly in the hole of the Whip Mix screw-type index. Tighten the finger screw on the front edge of the index (see Figure 15).


Figure 11


Figure 12

Figure 13

Figure 14


Figure 15

Whip Mix Quick Mount Index

Remove the mounting plate from the lower member. Place the Whip Mix Quick Mount index over the magnetic mounting base with the transfer assembly hole to the front of the articulator. With the numbers on clamps #1 and #2 in the upright position, secure the vertical shaft of the bitefork assembly in the hole of the Whip Mix Quick Mount index. Tighten the finger screw on the front edge of the index (see Figure 16).


Figure 16

Mounting The Maxillary Cast

Attach a mounting plate to the upper member. Be sure that the incisal pin is at the zero position and that the upper member is level and parallel to the table top.

Place the maxillary cast in the bite registration index on the bitefork and close the articulator to check for adequate clearance between the base of the cast and the mounting plate. Mount the cast with stone to the mounting plate (see Figure 17). Although it is not necessary, some clinicians prefer to use a maxillary cast support to prevent torque during the mounting procedure. Once the stone has set, remove the transfer jig and replace the incisal table in the articulator.

Figure 17

Mounting The Mandibular Cast

Attach a mounting plate to the lower member. A centric relation checkbite record is used to mount the mandibular cast. With the maxillary cast attached to the articulator and the centric latch engaged, invert the articulator and place the centric relation checkbite record between the maxillary and mandibular casts (see Figure 18). Stabilize the position of the casts with either a rubber band or sticky wax. Adjust the incisal pin to accommodate the increased vertical distance caused by the thickness of the centric relation bite record. Double check to be certain that the condyles are seated against the rear and tops of the fossae, and using stone, proceed to mount the mandibular cast to the mounting plate on the articulator (see Figure 19).

Figure 18


Figure 19

V. Slidematic Accessories

The Slidematic uses either a dentulous or edentulous bitefork. The Slidematic Facebow can be used with other articulators and indexes are available for the Hanau™ and Whip Mix® models.

Standard Biteforks

Indexes

VI. Reordering

Standard Biteforks

Part No. 200029 Dentulous
200030 Edentulous

Articulator Indexes

Part No. 200080 Denar®
200082 Hanau™
200088 Whip Mix®
200088Q Whip Mix® QuickMount

Other

Part No. 200009-1 Transfer Jig
Assembly
200007 Bitefork/Shaft
Assembly
101217 Maxillary Cast
Support

Denar® and Whip Mix® are registered trademarks and Hanau™ is a trademark of Whip Mix Corporation.

VII. Warranty

Whip Mix Corporation warrants the Slidematic Facebow to be free from defects in material and/or workmanship for a period of one year. In the event of a defect, please notify the factory in writing of the defect prior to returning the instrument. Whip Mix Corporation will, at its option, either repair, replace or issue credit for such defects.

Because Whip Mix Corporation is continually advancing the design of its products and manufacturing methods, it reserves the right to improve, modify or discontinue products at any time, or to change specifications or prices without notice and without incurring obligations.

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© 2012 Whip Mix Corporation
Denar® and logo are registered trademarks
Denar® y Whip Mix® son marcas comerciales registradas y Hanau™ es una marca comercial de Whip Mix Corporation.
Whip Mix Europe GmbH