

AEROETCHER™ BLASTER

Great for aluminium oxide air-abrasion of metal, resin and porcelain prior to adhesive procedures. Trigger action makes it equally accessible to right- and left-handed operators (no awkward button to push).

Like most intraoral blasters, the AeroEtcher's head unscrews for autoclaving. But if you've been blasting posterior restorations, sterilising just the head may not be enough to prevent cross infection.

So with the exception of its abrasive reservoir, the entire blaster is fully auto-clavable. (A quick disconnect at the back of the handpiece allows snap removal.) Operates on line pressure as low as 40 psi.



Blast inside crowns before cementing. Every cement (even zinc phosphate) provides more retention when you blast immediately before luting.

- AeroEtcher™ Intraoral Blaster
REF: D670
RRP **£134.00 + vat**
OFFER PRICE **£99.00 + vat**
- Aluminium Oxide Abrasive
REF: D675 **£22.60 + vat**
- Quick Connect fitting to plug into existing air outlet
REF: D671 **£23.50 + vat**



Blast before intraoral repairs. Triples the strength of composite-to-metal bonds. (The AeroEtcher isn't designed for caries removal or abrasion of tooth structure.)

waterpik NO 1 DENTISTS CHOICE DEEP CLEANING ORAL CARE

NEW Cordless, Handheld Rechargeable Irrigator

- Sleek, cordless one handed design
- Convenient
- Rechargeable, no batteries needed
- Clinically proven to reduce gingivitis
- Ideal for cleaning around orthodontic brackets
- Ideal for implant aftercare
- A great practice revenue builder
- Two tips for gentle cleaning
- Two tips for deep cleaning
- Tips rotate 360°
- Easy to fill reservoir may be used with just water or medicated solutions



MARKETING OPPORTUNITY
Why not give one to patients with a more expensive treatment plan. They'll think your the best!

RRP **£17.84 + vat**
REF: WP - 360



Prestige PRACTICALITIES

HEAD OFFICE 7 Oxford Place, Bradford, West Yorkshire BD3 0EF Tel: 01274 721567 Fax: 01274 304237

LONDON OFFICE 29 Chiltern Street, London W1U 7PL Tel: 020 7486 3881 Fax: 020 7224 2084

WEB www.prestige-dental.co.uk EMAIL info@prestige-dental.co.uk FREEPHONE 0800 591175 FREEFAX 0800 0282714

TAKE ADVANTAGE OF OUR EXCELLENT PRICES...

PLUS we will give you a further 10% off the first 10 items ordered from this brochure (see back page for details)

NEW

EARTHSTONE ORTHO... ORTHO-SIZED GYPSUM COMPOSITE

Provides twice as much stone per mix for larger orthodontic models in a bright white colour.

Simply shake and pour to create an instant silky smooth bubble free model.

New Ortho-Sized Earthstone
REF: 70-88013

RRP **£19.80 + vat** bx 12

The scoop, the dust, the mess, the waste, the clean-up, the long setting time - all are history. Discover the greatest innovation in model preparation since the spatula. Save valuable chair time and improve productivity for dentist, staff and patient. Many tasks that took multiple office visits can now be accomplished in one appointment.



PLUS NEW

ORTHO STUDY MODEL MOLDS

Saves you time... no model trimming, easy to clean. Molds have a built-in lingual contoured edge, with a gap for tray handles. 100% silicone construction for a long life.

Available in a set that contains 3 graduated upper and lower molds.

Six Ortho Study Model Molds in total

@ **£35.70 per set + vat**
REF: 70-82815



10% DISCOUNT IS AVAILABLE ON THE PURCHASE OF THE FIRST 10 PRODUCTS FROM THIS BROCHURE

QUANTITY MAY BE MADE UP OF INDIVIDUAL PRODUCTS OR MULTIPLES OF ONE PRODUCT. FOC PRODUCTS DO NOT COUNT AS 1 OF THE 10 PRODUCTS REQUIRED TO ACHIEVE DISCOUNT.

Please complete and return

PRODUCT	CODE	QUANTITY
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		

DENTIST NAME _____
ADDRESS _____
POSTCODE _____
TEL NO _____
EMAIL _____

Please return to
Prestige Dental, 7 Oxford Place, Bradford,
West Yorkshire BD3 0EF

Or Priority Fax Order to **01274 304237**



BRADFORD OFFICE: 7 OXFORD PLACE, BRADFORD, WEST YORKSHIRE BD3 0EF
TEL: 01274 721567 FAX: 01274 304237 EMAIL: info@prestige-dental.co.uk WEB: www.prestige-dental.co.uk

HOW ABOUT THIS... PURCHASE 5 BOXES OF ORTHO EARTHSTONE AND RECEIVE 1 SET OF ORTHO STUDY MOLDS FREE - PLUS YOUR POSSIBLE EXTRA 10% OFF THE ORTHO EARTHSTONE (see back page for details)

NEW LASCHAL INSTRUMENTS - WHEN ONLY THE BEST WILL DO!

LASCHAL BLEACHING TRAY AND RETAINER SCISSORS

Designed to easily and accurately trim vacuum formed trays and retainers etc. up to 0.02" thick.

REF: 70-86311
RRP **£17.75 + vat**



NEW LASCHAL INSTRUMENTS - WHEN ONLY THE BEST WILL DO!

LASCHAL SURGICAL BRACKET FORCEPS

Specifically designed for precision handling and delivery of orthodontic eruption appliances to surgically exposed cuspids.

REF: 70-86381
RRP **£133.35 + vat**



LASCHAL INSTRUMENTS... WHEN ONLY THE BEST WILL DO

FOR A FULL LASCHAL INSTRUMENTS CATALOGUE CALL TODAY **01274 721567**

ORDER FREEPHONE **0800 591175**

ORDER FREEFAX **0800 0282714**

EMAIL: info@prestige-dental.co.uk

INTRA ORAL PHOTOGRAPHY MIRRORS AT A REASONABLE PRICE!



Narrow Buccal REF: 70-61910
Wide Lingual REF: 70-61911
Large Occlusal REF: 70-61912
Child Occlusal REF: 70-61913

■ Rhodium coated - for extra long life ■ Autoclavable ■ Double-sided
JUST £28.50 + vat with your optional 10% extra off

LASCHAL - EH SUTURE SCISSORS

PROBLEM
 A problem exists when a patient returns with sutures buried deeply within edematous or hypertrophic (swollen) tissue. Locating and severing the buried suture with traditional surgical or Littauer scissors can cause additional trauma, pain bleeding and general discomfort to the patient.

SOLUTION
 EH Suture Scissors provide a significantly improved method of locating and severing buried sutures without traumatizing swollen, sensitive tissue. They eliminate the need to enter from beneath the suture or elevate the tag above the knot for access.

REF: 70-86312 RRP £53.99 + vat



LASCHAL - SPRING LOCK NEEDLE HOLDERS



Heavier gauge Spring-lock Needle Holders offer the weight and feel of a heat-treated instrument in a Laschal Featherlight design.

REF: 70-86315 RRP £140.65 + vat 5 1/2"
 REF: 70-86314 RRP £140.65 + vat 7 1/2"

DIAMOND INTERPROXIMAL TOOTH STRIPPER BY NTI

With the added advantage of a serrated edge

Diamond Finishing Strips, with a serrated edge to cut through any flash in interdental regions. Suitable for use on amalgam, composite, acrylics etc.

Available in three grit sizes and colour coded - Yellow - Superfine Red - Fine Blue - Medium

AND THEY ARE AUTOCLAVABLE!

AVAILABLE @ £29.90 pk/10 + vat



FLURO-TRAYS - FLUORIDE TRAYS

Anatomically contoured fluoride application trays packaged in 50 upper trays white and 50 lower trays blue. Available in 3 sizes.

Fluro Trays Small REF: 26200

Fluro Trays Medium REF: 26205

Fluro Trays Large REF: 26210


£14.45 pk/50 + vat



HURRAY!

STERILISATION POUCHES WITH AN INTERNAL PROCESSING INDICATOR BUILT-IN!


*ELIMINATE THE NEED FOR SEPARATE INDICATORS



31610
 (90MM x 135MM)
 Implant, Dam Clamp
 £16.75 + vat (200)



31613
 (55MM x 100MM)
 Burs
 £12.05 + vat (200)



31614
 (135MM x 165MM)
 Impression Tray, Bur Block
 £19.95 + vat (200)



31611
 (70MM x 230MM)
 Handpiece, Hand Instrument
 £13.95 + vat (200)



31612
 (90MM x 230MM)
 Handpiece, Hand Instrument
 £13.95 + vat (200)



31616
 (83MM x 305MM)
 Hand Instrument (Long)
 £17.50 + vat (200)



31617
 (108MM x 305MM)
 Surgical, Orthodontic
 £21.25 + vat (200)



31619
 (135MM x 305MM)
 Multiple Instrument
 £25.75 + vat (200)



31618
 (203MM x 305MM)
 Multiple Instrument, Small Cassette
 Stericage, AA and AB SteriSystem Cassette
 £35.20 + vat (200)



31615
 (216MM x 381MM)
 Multiple Instrument, AA SteriSystem Cassette
 £48.75 + vat (200)



31643
 (265MM x 405MM)
 Large Cassettes, BB SteriSystems
 £69.85 + vat (200)



31647
 (335MM x 510MM)
 Large Cassettes, (I.M.S. Type)
 £104.45 + vat (200)

PeelVue+ Pouches are easy to peel and seal. They feature a patented Internal-Processing-Indicator (IPI), and are printed with lead-free inks on medical-grade, virgin, non-recycled paper. After removing the clear film, the paper side of the pouch can function as a sterilized tray liner.



THE SMART POUCH Here is what PeelVue+™ Sterilization POUCH OFFER ...FEATURE A PATENTED INTERNAL PROCESSING INDICATOR...

Internal and External Indicators built directly into each pouch...

A solid seal on the peripheral edge of the pouch enhancing the strength of the pouch and preventing instruments from penetrating through.

A pre-applied sealing tape with a unique non-adhesive edge for easy peel and seal of the pouch - even with gloves. The sealing tape is designed into the pouch.

PeelVue Pouches are delivered in pre-cut lengths designed with Dentistry in mind. No racks, no sealing machines, no tape, no waste. Save space, time and money.

Vents located on either side of the internal indicator to allow the same ambient conditions inside and outside the pouch.

A Chevron Seal preventing instruments from touching the indicator and falsifying a proper reading.

Clear plastic for easy viewing of contents in the pouch - particularly important when viewing colour coded instruments and the internal indicator.



THE DENTAL ADVISOR



TAKE ADVANTAGE OF OUR EXCELLENT PRICES... PLUS YOUR EXTRA 10% OPTION (see back page for details)

PULPDENT WET BOND PIT & FISSURE SEALANT 'EMBRACE'




Light Cure, Fluoride Releasing

The first wet bonding pit and fissure sealant 'Embrace' seals against microleakage, is margin free so no chipping or staining. No drying or bonding agents are required.

Simply clean and etch, apply Embrace and Light Cure.

Bonds in the presence of moisture (even saliva).



Cleaned, etched and rinsed preparation. Leave tooth surface slightly moist.



EMBRACE natural shade applied to moist tooth and light cured. Also available in off-white.

Kit contains 4 x 1.2ml syringes and 20 applicator tips

Natural Shade REF: EMS

PRICE £41.95 + vat

Off-white Shade REF: EMSW

PRICE £45.95 + vat

ORDER FREEPHONE 0800 591175

ORDER FREEFAX 0800 0282714

EMAIL: info@prestige-dental.co.uk

PULPDENT ETCHRITE

A viscous 38% phosphoric acid. It stays where it is placed!

REF: ETCH

Kit contains 4 x 12ml syringes and 8 prebent tips

£14.30 + vat Kit

Alternatively available in 12ml syringes

REF: ET12 @£9.90 + vat

PULPDENT ETCHRITE PRE BENT TIPS

REF: ETL £18.80 + vat pk 100



PULPDENT PORCELAIN PREPARATION KIT



Prepare Porcelain for bonding at a realistic price.

REF: PPK

One economical kit containing Porcelain Etch Gel, Drying Agent and Silane.
 £26.60 + vat

WOULD YOU LIKE A FULL PULPDENT PRODUCT CATALOGUE? CALL 01274 721567

TRA-TEN IMPRESSION TRAYS


Ideal for alginate impressions

RRP £8.70 + vat pk 12

OFFER PRICE £6.90 + vat pk 12

Plus your extra 10% option

ORDER 20 PACK ANY ASSORTED SIZES AND WE WILL SEND YOU 1KG CADCO 'IMAGE' DUST FREE, FAST SET ALGINATE FREE OF CHARGE



HAVE YOU EVER CONSIDERED A CADCO ALGINATE MIXER? CALL 01274 721567 FOR A DEMO.