



Pasten, Emulsion

Das richtige Poliermittel ist entscheidend für eine erfolgreiche Politur. Um effektiv die unterschiedlichen Werkstoffe schrittweise zu polieren, werden bei den HATHO Polistar Polierpasten ausgewählte Wachse und hochwertige Mineralien miteinander kombiniert. Für eine optimale Zuordnung der Anwendungsgebiete werden den Polistar Rezepturen natürliche Farben beigemischt. Der hohe Füllgrad an Mineralien zeichnet das HATHO Polistar Sortiment aus. Als krönenden Abschluss einer erfolgreichen Hochglanzpolitur findet die Polistar Emulsion ihre Anwendung.



Compounds, Emulsion

Using the correct polishing compound is the most important thing in the world of polishing.

To polish different materials step by step quickly and effectively, HATHO offers its Polistar range. Selected waxes and high quality minerals are optimally combined.

Polistar products are distinguished by a very high fill density of the minerals. For the very last high shining step the liquid Polistar Emulsion is used.



Art.Nr. art.No.	ID.Nr. ref. No.	Anwendungen Applications
 ~300g / ~90-100g	Polistar Braun 300g Polistar Brown 300g Polistar Braun 100g Polistar Brown 100g 313177 313176	Vorpolieren von Edel- und Nichtedelmetalllegierungen Pre-polishing of precious and non-precious metal
	Polistar Lintygrau 300g Polistar Lintygrey 300g Polistar Lintygrau 90g Polistar Lintygrey 90g 313175 313168	Vorpolieren von Sterling Silber, Dentalkunststoffen und Komposite sowie Glanzmatieren von Keramik Pre-polishing of sterling silver, acrylic dentures, composites as well as shine matting of ceramic
	Polistar Weiss 300g Polistar White 300g Polistar Weiss 100g Polistar White 100g 313169 313162	Glanzpolitur von Platin, Nichtedelmetallen. CrCo-Modellguß, Thermoplaste Shine polishing of Platinum, non-precious metal, CrCo-framework, thermoplastic resin
	Polistar Grün 300g Polistar Green 300g Polistar Grün 100g Polistar Green 100g 313173 313166	Hochglanzpolitur von Edel- und Nichtedelmetalllegierungen, CrCoW-Modellguß High Shine polishing of precious and non precious metal alloys, CrCoW-framework
	Polistar Blau 300g Polistar Blue 300g Polistar Blau 100g Polistar Blue 100g 313172 313165	Hochglanzpolitur von Platin, Nichtedelmetalle, CrCoMo-Modellguß, Hybrid/Nano Komposite High Shine polishing of Platinum, non-precious metal, CrCoMo- framework, hybride/nano composites
	Polistar Rot 300g Polistar Red 300g Polistar Rot 100g Polistar Red 100g 313170 313163	Hochglanzpolitur von 18 Karat Gold, AuPt-Goldlegierungen, Glanzpolitur von Titan High Shine polishing of 18 carat gold, high noble gold alloys. Shine polishing of titanium
	Polistar Beige 300g Polistar Beige 300g Polistar Beige 100g Polistar Beige 100g 313174 313167	Hochglanzpolitur von Dentalkunststoffen und Komposite (Hartwaxschutzschicht) High Shine polishing of acrylic dentures and composites (Wax-Protection-Layer)
	Polistar Rosa 300g Polistar Rose 300g Polistar Rosa 100g Polistar Rose 100g 313171 313164	Hochglanzpolitur von 14 Karat Gold, Platin, Hybrid Komposite und Thermoplaste High Shine polishing of 14 carat gold, Platinum, hybride composites and thermoplastic resin
	Polistar Diamant 2,5g Polistar Diamond 2,5g Polistar Diamant 1,5g Polistar Diamond 1,5g 313079 313085	Hochglanzpolitur von Edel- und Nichtedelmetallen, CrCo-Modelguss, Titan, Nano Komposite und Dentalkeramik High Shine polishing of precious and non-precious metal, CrCo-framework, titanium, nano composites and dental ceramics
	Polistar Emulsion 125ml Polistar Emulsion 125ml Polistar Emulsion 50ml Polistar Emulsion 50ml 313139 313148	Finale Hochglanzpolitur von Edel- und Nichtedelmetallen, CrCo-Modelguss, Titan, Komposite und Dentalkeramik Final High Shine polishing of precious and non-precious metal, CrCo-framework, titanium, composites and dental ceramics