



Technically first class

RESTORATIVE ESSENTIALS

ETCHING AND CLEANING



PULPDENT ETCH-RITE

Etch-Rite is a soft, thixotropic blue gel with handling characteristics most preferred by clinicians.

It dispenses through small gage needles, stays where placed, washes off with ease and provides the optimal etch pattern on dentine and enamel surfaces to ensure mechanical retention of bonding agents, restorative resins and resin cements.

4 x 1.2ml syringes + 8 tips

Order ref: ETCH

£19.20



AIRSONIC MINI SANDBLASTER

Handy sandblaster for roughening of bonding surfaces of crowns, bridges, brackets, posts etc.

Order ref: 401080

£225.00

Adapters available for: Kavo, Sirona, W&H, BienAir, NSK, EMS, Lab (sold separately).



BISCO PORCELAIN ETCHANT

Porcelain Etchant improves the adhesive bond between porcelain and cement by increasing the surface area of the porcelain to be bonded. The effect is similar to sandblasting but is more controlled and effective. Porcelain Etchant will microscopically roughen the feldspathic porcelain or glass surface.

1 x 5g syringe + 25 tips
4% Hydrofluoric acid gel

Order ref: E-5702EP

£21.40

1 x 5g syringe + 25 tips
9.5% Hydrofluoric acid gel

Order ref: E-5707P

£21.40

BISCO UNI-ETCH

Uni-Etch w/BAC is a 32% phosphoric acid semi-gel etchant that is effective at quickly etching enamel and dentin. Uni-Etch w/BAC is polymer-thickened and effectively removes the smear layer, etches and demineralizes enamel and dentin. It is easily manipulated into small areas such as endodontically treated root canals and undercuts. This makes it more versatile than some thicker gel etchants for etching all restorative surfaces.

4 x 5g syringes + 60 tips

Order ref: E5502EBM

£41.45

1 x 30g bottle

Order ref: E5637EB

£32.10

1 x 30ml syringe

Order ref: E-56621P

£35.40



BONDING/PULP CAPPING

PULPDENT DENTASTIC UNO/DUO

Dentastic Uno is an effective, single component light cure adhesive for bonding to Dentine, enamel, porcelain, metal, composite and other resins. Uno is used for all direct bonding light cure applications Dentastic Duo is the dual cure catalyst for Dentastic Uno. Use Uno+Duo whenever self-cure or dual cure capability is indicated.

2 x 6ml Uno, 5ml Etch-Rite
+ 20 tips and 50 brush tips

Order ref: UNO

£155.00

3ml Bottle Duo Catalyst

Order ref: DUO

£70.00

6ml Uno, 3ml Duo, 5ml Etch-Rite
+ 20 tips and 50 brush tips

Order ref: UNDO

£145.00



BISCO ALL-BOND UNIVERSAL

As a universal adhesive it can be used with direct and indirect restorations and is formulated to be compatible with light-, dual- and self-cured materials. The versatility of All-Bond Universal makes it an indispensable part of any dental practice.

All Bond Universal Standard Kit
6ml All-Bond U, 5g syringe Select HV Etch

Order ref: B-72020K

£95.45

BISCO THERACAL LC

TeraCal LC is a light-cured, resin-modified calcium silicate. Its unique apatite stimulating ability makes it ideal for direct and indirect pulp capping and as a protective liner.

1 x 1g syringe + accessories

Order ref: H-3301P

£27.40

4 x 1g syringes + accessories

Order ref: H-33014P

£86.20



PULPDENT MIXING WELLS

These convenient disposable mixing wells are compatible with all solvents, bonding agents and dental materials. Available in 2-well and 4-well configurations.

Box of 480 2-well
configuration mixing wells

Order ref: MW-2

£57.45

Box of 420 4-well
configuration mixing wells

Order ref: MW-4

£57.45



LINING AND COMPOSITE

PULPDENT LIME-LITE ENHANCED

Lime-Lite™ Enhanced is an improved moisture-friendly formula specially formulated for use with adhesives, composites and conventional restorative materials. (Not a pulp cap)

4 x 1.2ml syringes + 8 tips
Order ref: **LIME**
£56.00

1 x 3ml syringe
Order ref: **LIME3**
£33.45



PULPDENT ACTIVA BIOACTIVE RESTORATIVE

ACTIVA dual cure BioACTIVE-RESTORATIVE delivers all the advantages of glass ionomers in a strong, resilient resin matrix that will not chip or crumble. It chemically bonds to teeth, seals against bacterial microleakage, releases more fluoride, is more bioactive than glass ionomers, and is more durable and fracture resistant than composites.

ACTIVA is the first bioactive composite with an ionic resin matrix, a shock-absorbing resin component and bioactive fillers that mimic the physical and chemical properties of natural teeth. It releases and recharges with calcium, phosphate and fluoride ions, delivering long term benefits for your patients.

+ BPA/Chemical Free



PULPDENT ACTIVA BIOACTIVE BASE LINER

ACTIVA dual cure BioACTIVE-BASE/LINER delivers all the advantages of glass ionomers in a strong, resilient resin matrix that will not chip or crumble.

1 x 5ml syringe + 20 tips
Order ref: **VB1**
£119.00

2 x 5ml Syringes + 20 tips
Order ref: **VB2**
£198.45



Starter Kits		Order ref:	Price:
A1	1 x 5ml syringe, dispensing gun + 20 tips	VRA1	£164.95
A2	1 x 5ml syringe, dispensing gun + 20 tips	VRA2	£164.95
A3	1 x 5ml syringe, dispensing gun + 20 tips	VRA3	£164.95
Single Refills			
A1	1 x 5ml syringe + 20 tips	VR1A1	£99.50
A2	1 x 5ml syringe + 20 tips	VR1A2	£99.50
A3	1 x 5ml syringe + 20 tips	VR1A3	£99.50
A3.5	1 x 5ml syringe + 20 tips	VR1A35	£99.50
Value Refills			
A1	2 x 5ml syringes + 40 tips	VR2A1	£179.50
A2	2 x 5ml syringes + 40 tips	VR2A2	£179.50
A3	2 x 5ml syringes + 40 tips	VR2A3	£179.50
A3.5	2 x 5ml syringes + 40 tips	VR2A35	£179.50

FINISHING & POLISHING

BAUSCH ARTICULATING PAPER

Bausch 40µ micro-thin articulating papers are extremely thin and tear resistant and are coated with liquid colours on both sides. These papers mark precisely because of the thin material; false or smear contacts can thus be avoided. The special coating with liquid colours facilitates accurate marking of all occlusal contacts or occlusal interferences. Moist occlusal surfaces such as gold, ceramic, polished metal or acrylic, which are difficult to examine, do not pose any problem. The special colour coating with liquid colours consists of many colour-filled microcapsules. Even the slightest masticatory pressure can cause the capsules to burst and thus release the distinctly visible colour. Also repeated marking is possible because of colour regeneration.

**Box of 200 Blue strips
Straight**

Order ref: BK09

£13.60

**Box of 200 Red strips
Straight**

Order ref: BK10

£13.60

**Box of 100 Blue sheets
100x70mm**

Order ref: BK11

£24.05

**Box of 100 Red sheets
100x70mm**

Order ref: BK12

£24.05



NTI PROXYSTRIPS

The ProxyStrip is especially suitable for the proximal design of fillings of all kinds as well as stripping, slicing and interproximal reduction. Particularly in the posterior region, the ProxyStrip is characterised by its simple and extremely space saving use. The angle adjustable heads make work much easier – even with limited opening.

**ProxyStrip Assortment –
Includes handle and 5 assorted
strips**

Order ref: FS5-KP-SO

£39.75

Refill pack of 10 Medium Strips

Order ref: FS5-MKP

£46.85

Refill pack of 10 Fine Strips

Order ref: FS5-FKP

£46.85

Refill pack of 10 Superfine Strips

Order ref: FS5-SFKP

£46.85

Refill pack of 5 Sawblade Strips

Order ref: FS5-XKP

£23.45

ProxyStrip Handle

Order ref: FSG-1

£21.00

MILLER FORCEPS

Stainless steel straight forceps designed for holding papers and foils during use.

1 Pair Miller Forceps

Order ref: BK132

£13.50



FINISHING & POLISHING

NTI SERRATED DIAMOND FINISHING STRIPS

These serrated diamond finishing strips have a unique saw edge and finishing strip combining two instruments in one.

Pack of 10 Medium grit strips
Order ref: FS3-MX

£43.80

Pack of 10 Fine grit strips
Order ref: FS3-FX

£43.80

Pack of 10 Superfine strips
Order ref: FS3-SFX

£43.80

Pack of 10 assorted strips
Order ref: FS3-X-SO

£43.80

NTI SERRATED AND PERFORATED FINISHING STRIPS

The new perforated and serrated diamond finishing strips adapt particularly flexibly to the surface of the tooth.

Pack of 10 Medium grit strips
Order ref: FX3-MXP

£43.80

Pack of 10 Fine grit strips
Order ref: FS3-FXP

£43.80

Pack of 10 Superfine grit strips
Order ref: FS3-SFXP

£43.80

Pack of 10 assorted strips
Order ref: FS3-XP-SO

£43.80



NTI FLEXIPOL COMPOSITE POLISHING DISCS

Highly flexible polishing discs in two diameters (8mm and 12mm) and four grit sizes for producing a natural lustre on all nanohybrid and standard composites.

25 of each disc + 4 Mandrels
Order ref: PC1295-SET

£53.50

Refill packs of 25 discs available for each size/grit

£21.45

PULPDENT SPARKLE DIAMOND POLISHING PASTE

Produces a glass-like high lustre finish on porcelain. Gold, composite and metal.

4 x 1.2ml Syringes
Order ref: SPARK

£36.95

1 x 3ml Syringe
Order ref: SPARK3

£22.00



NTI CERAGLAZE CERAMIC/ZIRCONIA POLISHERS

Ceraglaze polishers are manufactured with a special rubber binder which absorbs the heat from the surface during the polishing process. This avoids any overheating. Set includes three types. Green- Pre-polish, Blue- Smoothing, Yellow- High Shine in four shapes – Large Bullet, Small Bullet, Inverted Cone and Knife Edge Lens.

Set of 12 Polishers

Order ref: SET1588

£97.15


All polishers available in refill packs of 3 or 6 polishers.

Refill pack of 3

£31.10

Refill pack of 6

£61.75



NTI UNIVERSAL COMPOSITE SET

Includes all necessary rotary instruments for caries elimination, removal of old fillings, cavity preparation classes I, II, III, IV, V, contouring, finishing and polishing of composite fillings.

Set contains: 4 x Arkansas Stones, 7 x Diamond burs, 6 x Tungsten carbide burs, 3 x polishers.

Set of 20 instruments

Order ref: SET-1672

£63.65

All instruments available in individual packs.

NTI DIAGLOSS POLISHERS

Two step diamond treatment system for composites. Specially selected diamond particles binder allow for contouring and high gloss polishing of composites. The special rubber binder has the ability to absorb heat during polishing. This prevents damage at thin junctions between the tooth and the composite. Diagloss are also suitable for polishing natural teeth following bracket removal due to the ability to absorb heat. Set includes two types – Coarse grit for pre-polishing and High Shine in six shapes – Large Bullet, Small Bullet, Large Cup, Inverted Cone, Disc and Pear.


Set of 12 Polishers

Order ref: SET-1729

£51.25

All polishers available in refill pack of 6 polishers

£25.15



Prestige Dental Ltd 7 Oxford Place, Bradford, West Yorkshire, BD3 0EF

Tel: 01274 721567

Email: info@prestige-dental.co.uk • www.prestige-dental.co.uk

If you prefer not to receive our special offers please contact us to unsubscribe

All prices exclude VAT and postage and are subject to change without prior notification.

