

5 REASONS TO TRY TheraCem[®] Ca

Dual-Cured, Self-Adhesive Calcium Releasing Cement

Rx Only

1 Contains Calcium

Continuous calcium ion release.¹

2 Contains MDP

Strong bond to zirconia, metal and alumina substrates without primer.

3 Easy Clean-up

Gel phase after tack cure facilitates removal of excess cement.

4 Alkaline pH

Generates an alkaline pH in minutes², which promotes pulp vitality.³

5 Radiopaque

Visible on radiograph to easily distinguish from caries.



1. Gleave CM, Chen L, Suh BI. Calcium & fluoride recharge of resin cements. Dent Mater. 2016 (32S):e26.
2. New Self-adhesive Resin Cement With Alkaline pH. Chen L, Gleave C, Suh B, J Dent Res96(A):#286, 2017
3. T. Okabe, M. Sakamoto, H. Takeuchi, K. Matsushima. Effects of pH on Mineralization Ability of Human Dental Pulp Cells. Journal of Endodontics. Volume 32, Number 3, March 2006

MC-50211TB

For more information visit www.bisco.com

